



BOUTIQUE  
WHITENING

*Ask your dentist about...*

BOUTIQUE  
WHITENING

## WHITENING KIT BY DAY

6% Hydrogen Peroxide Gel - With Potassium Nitrate

Boutique - By Day is a 6% hydrogen peroxide gel, containing potassium nitrate which is a proven desensitiser. Worn for only one hour per day, this gel is designed for patients who want a fast, comfortable and convenient solution to whiter teeth.

This gel releases large quantities of oxygen free radicals (which whiten the teeth) over a short period of time, therefore whitening the teeth, both quickly and safely.

### Key features:

- Optimum performance - Highest percentage formula available
- Fast results, which are visible in less than a week
- PH Neutral • Viscous formula to eliminate wash out from saliva
- In-built desensitiser
- Water integrated gel to reduce dehydration
- Minimum wear time
- Amazing presentation
- Custom packaging available to the practice



## CUSTOM TRAYS

When using a quality whitening gel, it would be wasted without a premium quality custom whitening tray. We manufacture our own trays to ensure we are providing the best fit possible.

Our tray material is rigid enough that it will hold the gel in place, but flexible enough that it offers comfort for the patient. All of our margins are carved on the model, which is clearly visible in the tray. This keeps a tight seal between the tooth and the gingiva, therefore reducing irritation to the soft tissues.

The research shows that using a reservoir has little effect on the level of whitening achieved. However, with our viscous whitening gel, we offer 'the best of both'. We incorporate a small reservoir, which allows space for the gel, but also a custom fit to the rest of the tooth surface.

### Key features:

- Carved Margins and tight seal trays
- Rigid but flexible tray material
- Scalloped edges for patient comfort
- Small reservoir, but allowing a custom fit tray



## WHITENING KIT BY NIGHT

16% Carbamide Peroxide Gel - With Potassium Nitrate

By Night is a 16% carbamide peroxide gel, containing potassium nitrate, which is proven desensitiser. Worn for 4 hours, or overnight, this gel is designed for patients who would prefer to not have their custom trays in during the day and would prefer to whiten whilst they are asleep.

This gel releases oxygen free radicals (which whiten the teeth) over a longer period of time, due to its sustained release. This gradual approach will still give fantastic results, due the length of time the gel stays active.

### Key features:

- Optimum performance - Highest percentage formula available
- Fast results, which are visible in less than a week
- PH Neutral
- Viscous formula to eliminate wash out from saliva
- In-built desensitiser
- Water integrated gel to reduce dehydration
- Sustained release formula
- Amazing presentation
- Custom packaging available to the practice



Before



After

Ask your dentist about...

**BOUTIQUE**  
WHITENING